

Product datasheet for **GA205446**

Mouse Tmem38b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem38b CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 38B
Symbol:	Tmem38b
Components:	GA205446G1, Tmem38b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205446G2, Tmem38b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205446G3, Tmem38b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9DAV9</u>
Synonyms:	1600017F22Rik; AA437809; AV051057; D4Ert89e; mg33b; TRIC-B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)