

Product datasheet for **GA205447**

Mouse Pvr activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Pvr CRISPRa kit - CRISPR gene activation of mouse poliovirus receptor |
| Symbol: | Pvr |
| Synonyms: | 3830421F03Rik; CD155; D7Ert458e; HVED; mE4; necl-5; PVS; Taa1; Tage4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)