

Product datasheet for GA205468

Mouse Macir activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	D1Ert622e CRISPRa kit – CRISPR gene activation of mouse DNA segment, Chr 1, ERATO Doi 622, expressed
Symbol:	Macir
Components:	GA205468G1, Macir gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA205468G2, Macir gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA205468G3, Macir gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q8VEB3</u>
Synonyms:	AI987691; AW047450
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

