

Product datasheet for **GA205507**

Mouse Cd300lg activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cd300lg CRISPRa kit - CRISPR gene activation of mouse CD300 molecule like family member G
Symbol:	Cd300lg
UniProt ID:	<u>Q1ERP8</u>
Synonyms:	2310016B05Rik; Clm9; D11Ert736e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)