

Product datasheet for **GA205509**

Mouse Setd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Setd3 CRISPRa kit - CRISPR gene activation of mouse SET domain containing 3
Symbol:	Setd3
Components:	GA205509G1, Setd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205509G2, Setd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205509G3, Setd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q91WC0</u>
Synonyms:	2610102I01Rik; 2610305M23Rik; 6030490B17Rik; D12Ert771e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)