

Product datasheet for GA205513

Mouse Zfp410 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp410 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 410
Symbol:	Zfp410
Components:	GA205513G1, Zfp410 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205513G2, Zfp410 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205513G3, Zfp410 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BKX7</u>
Synonyms:	Apal; D12Ertd748e; Znf410
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

