

Product datasheet for **GA205523**

Mouse Dnlz activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dnlz CRISPRa kit - CRISPR gene activation of mouse DNL-type zinc finger
Symbol:	Dnlz
Components:	GA205523G1, Dnlz gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205523G2, Dnlz gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205523G3, Dnlz gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D113</u>
Synonyms:	1110034G11Rik; D2Bwg1335e; HEP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)