

Product datasheet for **GA205539**

Mouse Dlg3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dlg3 CRISPRa kit - CRISPR gene activation of mouse discs large MAGUK scaffold protein 3
Symbol:	Dlg3
Components:	GA205539G1, Dlg3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205539G2, Dlg3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205539G3, Dlg3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P70175</u>
Synonyms:	Dlgh3; mKIAA1232; SAP102
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)