

Product datasheet for **GA205622**

Mouse Phgr1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Phgr1 CRISPRa kit - CRISPR gene activation of mouse proline/histidine/glycine-rich 1
Symbol:	Phgr1
Components:	GA205622G1, Phgr1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205622G2, Phgr1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205622G3, Phgr1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8K0G7</u>
Synonyms:	1810007E14Rik; AA238765
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)