

Product datasheet for GA205627

Mouse Cyp3a41a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp3a41a CRISPRa kit – CRISPR gene activation of mouse cytochrome P450, family 3, subfamily a, polypeptide 41A
Symbol:	Cyp3a41a
Components:	GA205627G1, Cyp3a41a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205627G2, Cyp3a41a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205627G3, Cyp3a41a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9JMA7</u>
Synonyms:	Cyp3a41
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

