

Product datasheet for **GA205649**

Mouse **Dnal4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Dnal4 CRISPRa kit - CRISPR gene activation of mouse dynein, axonemal, light chain 4
Symbol:	Dnal4
UniProt ID:	<u>Q9DCM4</u>
Synonyms:	D15Ert424e; Dnalc4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)