

Product datasheet for **GA205661**

Mouse Pabpn1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pabpn1 CRISPRa kit - CRISPR gene activation of mouse poly(A) binding protein, nuclear 1
Symbol:	Pabpn1
Components:	GA205661G1, Pabpn1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205661G2, Pabpn1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205661G3, Pabpn1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8CCS6
Synonyms:	mPABII; PAB2; Pabp3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)