

Product datasheet for **GA205695**

Mouse H1f9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hils1 CRISPRa kit - CRISPR gene activation of mouse histone H1-like protein in spermatids 1
Symbol:	H1f9
Components:	GA205695G1, H1f9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205695G2, H1f9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205695G3, H1f9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9QYL0</u>
Synonyms:	H1.9; TISP64
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)