

Product datasheet for **GA205726**

Mouse Nap1I3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Nap1I3 CRISPRa kit - CRISPR gene activation of mouse nucleosome assembly protein 1-like 3
Symbol:	Nap1I3
Components:	GA205726G1, Nap1I3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205726G2, Nap1I3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205726G3, Nap1I3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q794H2
Synonyms:	MB20
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)