

Product datasheet for **GA205780**

Mouse Ap3m1 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Ap3m1 CRISPRa kit - CRISPR gene activation of mouse adaptor-related protein complex 3, mu 1 subunit |
| Symbol: | Ap3m1 |
| UniProt ID: | Q9JKC8 |
| Synonyms: | 1200013D09Rik; C78982; R75378 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)