

Product datasheet for **GA205793**

Mouse Paxip1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Paxip1 CRISPRa kit - CRISPR gene activation of mouse PAX interacting (with transcription-activation domain) protein 1
Symbol:	Paxip1
Components:	GA205793G1, Paxip1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205793G2, Paxip1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205793G3, Paxip1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6NZQ4</u>
Synonyms:	D5Ert149e; PTIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)