

Product datasheet for **GA205794**

Mouse **Pdzrn3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Pdzrn3 CRISPRa kit - CRISPR gene activation of mouse PDZ domain containing RING finger 3
Symbol:	Pdzrn3
UniProt ID:	Q69ZS0
Synonyms:	1110020C07Rik; AI429718; AL023082; LNX3; SEMACAP3; Semcap3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)