

Product datasheet for **GA205804**

Mouse Alyref2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Alyref2 CRISPRa kit – CRISPR gene activation of mouse Aly/REF export factor 2
Symbol:	Alyref2
Components:	GA205804G1, Alyref2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205804G2, Alyref2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205804G3, Alyref2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9JJW6
Synonyms:	C130042O11Rik; REF2; Refbp2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

