

Product datasheet for **GA205833**

Mouse Psmg1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Psmg1 CRISPRa kit - CRISPR gene activation of mouse proteasome (prosome, macropain) assembly chaperone 1
Symbol:	Psmg1
UniProt ID:	Q9JK23
Synonyms:	AW552102; Dscr2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)