

Product datasheet for **GA205853**

Mouse Ptbp2 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Ptbp2 CRISPRa kit - CRISPR gene activation of mouse polypyrimidine tract binding protein 2 |
| Symbol: | Ptbp2 |
| Components: | GA205853G1, Ptbp2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205853G2, Ptbp2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205853G3, Ptbp2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>Q91Z31</u> |
| Synonyms: | brPTB; nPTB; Ptb2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)