

Product datasheet for **GA205867**

Mouse Pals1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mpp5 CRISPRa kit - CRISPR gene activation of mouse membrane protein, palmitoylated 5 (MAGUK p55 subfamily member 5)
Symbol:	Pals1
Components:	GA205867G1, Pals1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205867G2, Pals1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205867G3, Pals1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9JLB2
Synonyms:	3830420B02Rik; AI255216; AI644496; Pals1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)