

Product datasheet for **GA205873**

Mouse Tspan5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tspan5 CRISPRa kit - CRISPR gene activation of mouse tetraspanin 5
Symbol:	Tspan5
UniProt ID:	P62080
Synonyms:	2810455A09Rik; 4930505M03Rik; NET-4; Tm4sf9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)