

Product datasheet for **GA205903**

Mouse Zfp113 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp113 CRISPRa kit – CRISPR gene activation of mouse zinc finger protein 113
Symbol:	Zfp113
Components:	GA205903G1, Zfp113 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205903G2, Zfp113 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205903G3, Zfp113 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	4732456B05Rik; mKIAA4229
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

