

## Product datasheet for **GA205940**

### Mouse **Cyp3a25** activation kit by **CRISPRa**

#### Product data:

<b>Product Type:</b>	CRISPRa Kits
<b>Description:</b>	Cyp3a25 CRISPRa kit - CRISPR gene activation of mouse cytochrome P450, family 3, subfamily a, polypeptide 25
<b>Symbol:</b>	Cyp3a25
<b>Components:</b>	GA205940G1, Cyp3a25 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205940G2, Cyp3a25 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205940G3, Cyp3a25 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
<b>UniProt ID:</b>	<a href="#">O09158</a>
<b>Synonyms:</b>	AI327008
<b>Format:</b>	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

