

Product datasheet for **GA205951**

Mouse *Trpc4ap* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trpc4ap CRISPRa kit - CRISPR gene activation of mouse transient receptor potential cation channel, subfamily C, member 4 associated protein
Symbol:	Trpc4ap
UniProt ID:	Q9JLV2
Synonyms:	4833429F06Rik; D2Erttd113e; mFLJ00177; Trrp4ap; truss
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)