

## Product datasheet for **GA205968**

### Mouse Tspan3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tspan3 CRISPRa kit - CRISPR gene activation of mouse tetraspanin 3
Symbol:	Tspan3
UniProt ID:	<u><a href="#">Q9QY33</a></u>
Synonyms:	1700055K04Rik; TM4-A; Tm4sf8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)