

Product datasheet for GA205977

Mouse Cyp2d22 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp2d22 CRISPRa kit - CRISPR gene activation of mouse cytochrome P450, family 2, subfamily d, polypeptide 22
Symbol:	Cyp2d22
Components:	GA205977G1, Cyp2d22 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA205977G2, Cyp2d22 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA205977G3, Cyp2d22 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
Synonyms:	2D22
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

