

Product datasheet for **GA205983**

Mouse Dynll1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dynll1 CRISPRa kit - CRISPR gene activation of mouse dynein light chain LC8-type 1
Symbol:	Dynll1
UniProt ID:	P63168
Synonyms:	Dlc8; Dncl1; Pin
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)