

Product datasheet for **GA205988**

Mouse Pkp3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pkp3 CRISPRa kit - CRISPR gene activation of mouse plakophilin 3
Symbol:	Pkp3
UniProt ID:	<u>Q9QY23</u>
Synonyms:	2310056L12Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)