

Product datasheet for **GA205989**

Mouse Kcnip3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Kcnip3 CRISPRa kit - CRISPR gene activation of mouse Kv channel interacting protein 3, calsenilin |
| Symbol: | Kcnip3 |
| UniProt ID: | <u>Q9QXT8</u> |
| Synonyms: | 4933407H12Rik; AI413860; Csen; DREAM; KChIP3; R74849 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)