

Product datasheet for **GA206004**

Mouse **Get3** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Asna1 CRISPRa kit - CRISPR gene activation of mouse arsa arsenite transporter, ATP-binding, homolog 1 (bacterial)
Symbol:	Get3
UniProt ID:	Q54984
Synonyms:	1810048H22Rik; ArsA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)