

Product datasheet for **GA206024**

Mouse Ripk3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ripk3 CRISPRa kit - CRISPR gene activation of mouse receptor-interacting serine-threonine kinase 3
Symbol:	Ripk3
UniProt ID:	<u>Q9QZL0</u>
Synonyms:	2610528K09Rik; AW107945; Rip3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)