

Product datasheet for **GA206031**

Mouse **Kcnd3** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Kcnd3 CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, Shal-related family, member 3 |
| Symbol: | Kcnd3 |
| UniProt ID: | Q9Z0V1 |
| Synonyms: | AW045978; Kcnd3; Kv4.3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)