


Product datasheet for **GA206039**

Mouse Clec4n activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Clec4n CRISPRa kit - CRISPR gene activation of mouse C-type lectin domain family 4, member n
Symbol:	Clec4n
Components:	GA206039G1, Clec4n gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206039G2, Clec4n gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206039G3, Clec4n gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9JKE4
Synonyms:	Clec6a; Clec4n; Nkcl
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)