

Product datasheet for **GA206047**

Mouse Klk4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Klk4 CRISPRa kit - CRISPR gene activation of mouse kallikrein related-peptidase 4 (prostase, enamel matrix, prostate)
Symbol:	Klk4
Components:	GA206047G1, Klk4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206047G2, Klk4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206047G3, Klk4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Z0M1</u>
Synonyms:	ESMP1; KLK-L1; Prss17; PSTS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)