

Product datasheet for GA206055

Mouse Phax activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Phax CRISPRa kit - CRISPR gene activation of mouse phosphorylated adaptor for RNA export
Symbol:	Phax
Components:	GA206055G1, Phax gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206055G2, Phax gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206055G3, Phax gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9JJT9</u>
Synonyms:	2810055C14Rik; 4933427L19Rik; AU018701; AU018854; D18Ert65e; p55; Rnuxa
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

