

Product datasheet for **GA206107**

Mouse Trpm5 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Trpm5 CRISPRa kit - CRISPR gene activation of mouse transient receptor potential cation channel, subfamily M, member 5 |
| Symbol: | Trpm5 |
| UniProt ID: | Q9JH7 |
| Synonyms: | 9430099A16Rik; LTrpC-5; Ltrpc5; Mtr1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)