

## Product datasheet for **GA206151**

### Mouse Lenep activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Lenep CRISPRa kit - CRISPR gene activation of mouse lens epithelial protein
Symbol:	Lenep
UniProt ID:	<a href="#">Q9WVB6</a>
Synonyms:	AU018355; Lep503
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)