

Product datasheet for **GA206157**

Mouse *Icmt* activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | <i>Icmt</i> CRISPRa kit - CRISPR gene activation of mouse isoprenylcysteine carboxyl methyltransferase |
| Symbol: | <i>Icmt</i> |
| UniProt ID: | Q9EQK7 |
| Synonyms: | 1700008E11Rik; C80758; Gm13095; HSTE14; OTTMUSG00000010406; pcCMT; PPMT; STE14 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)