

Product datasheet for **GA206159**

Mouse Mrps31 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mrps31 CRISPRa kit - CRISPR gene activation of mouse mitochondrial ribosomal protein S31
Symbol:	Mrps31
Components:	GA206159G1, Mrps31 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206159G2, Mrps31 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206159G3, Mrps31 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q61733</u>
Synonyms:	1500002D03Rik; AA960481; Imogen38; Imogen44
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)