

## Product datasheet for **GA206165**

### Mouse Terf2ip activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Terf2ip CRISPRa kit - CRISPR gene activation of mouse telomeric repeat binding factor 2, interacting protein
Symbol:	Terf2ip
UniProt ID:	<a href="#">Q91VL8</a>
Synonyms:	Rap1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)