


Product datasheet for **GA206179**

Mouse P2ry4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2ry4 CRISPRa kit - CRISPR gene activation of mouse pyrimidinergetic receptor P2Y, G-protein coupled, 4
Symbol:	P2ry4
Components:	GA206179G1, P2ry4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206179G2, P2ry4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206179G3, P2ry4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9JJS7</u>
Synonyms:	P2Y4; P2Y4R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)