

Product datasheet for **GA206181**

Mouse Cltrn activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cltrn CRISPRa kit - CRISPR gene activation of mouse collectrin, amino acid transport regulator
Symbol:	Cltrn
Components:	GA206181G1, Cltrn gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206181G2, Cltrn gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206181G3, Cltrn gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9ESG4</u>
Synonyms:	0610008J07Rik; NX-17; NX17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)