

Product datasheet for **GA206190**

Mouse Tmem183a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem183a CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 183A
Symbol:	Tmem183a
Components:	GA206190G1, Tmem183a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206190G2, Tmem183a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206190G3, Tmem183a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9JJB9
Synonyms:	1300007B12Rik; 1700020N19Rik; Tmem183
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)