

Product datasheet for **GA206193**

Mouse Kcne3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kcne3 CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, Isk-related subfamily, gene 3
Symbol:	Kcne3
UniProt ID:	Q9WTW2
Synonyms:	2210017H05Rik; MiRP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)