

Product datasheet for **GA206203**

Mouse Piwil2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Piwil2 CRISPRa kit - CRISPR gene activation of mouse piwi-like RNA-mediated gene silencing 2
Symbol:	Piwil2
UniProt ID:	<u>Q8CDG1</u>
Synonyms:	mili; Piwil11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)