


## Product datasheet for **GA206229**

### Mouse Kcne4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Kcne4 CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, Isk-related subfamily, gene 4
Symbol:	Kcne4
UniProt ID:	<a href="#">Q9WTW3</a>
Synonyms:	AW487937; MiRP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)