

Product datasheet for **GA206239**

Mouse Zfp318 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp318 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 318
Symbol:	Zfp318
Components:	GA206239G1, Zfp318 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206239G2, Zfp318 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206239G3, Zfp318 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q99PP2
Synonyms:	2610034E08Rik; D530032D06Rik; TZF; Znf318
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

