

Product datasheet for **GA206258**

Mouse *Sh3kbp1* activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Sh3kbp1 CRISPRa kit - CRISPR gene activation of mouse SH3-domain kinase binding protein 1 |
| Symbol: | Sh3kbp1 |
| UniProt ID: | <u>Q8R550</u> |
| Synonyms: | 1200007H22Rik; 1700125L08Rik; 5830464D22Rik; AI447724; Cin85; IN85; Ruk; Seta |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)