

Product datasheet for **GA206268**

Mouse *Srrm3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Srrm3</i> CRISPRa kit - CRISPR gene activation of mouse serine/arginine repetitive matrix 3
Symbol:	<i>Srrm3</i>
UniProt ID:	Q80WV7
Synonyms:	2900083111Rik; AB041661; SRm300-like; <i>Srrm2l</i>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)